

Product Datasheet

NCBP2 antibody (orb693076)

Description	Rabbit polyclonal antibody to NCBP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NCBP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18001
Uniprot ID	P52298
NCBI	NP_001036005.1, NP_031388.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.