



Product Datasheet

KLHL1 Antibody (orb69322)

Description	Mouse monoclonal antibody to KLHL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of human KLHL1 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	82kDa
Entrez	57626
Dilution Range	ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 C...



Western blot analysis of human KLHL1 rec...