



Product Datasheet

TNFRSF11B Antibody (orb69337)


Description	Mouse monoclonal antibody to TNFRSF11B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC
Immunogen	Purified recombinant fragment of human TNFRSF11B expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	60kDa
Entrez	4982
Dilution Range	ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HL-60 cells ...



Western blot analysis of HEK293 (Lane1) ...



Immunofluorescence analysis of HL-60 cel...