



Product Datasheet

LHX2 Antibody (orb69405)


Description	Mouse monoclonal antibody to LHX2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of human LHX2 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	44kDa
Entrez	9355
Dilution Range	ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of human LHX2 (AA:...



Immunofluorescence analysis of HeLa cell...



Line graph illustrates about the Ag-Ab r...