


Product Datasheet

SHIP1 (phospho-Tyr1020) antibody (orb6942)

Description	SHIP1 (phospho-Tyr1020) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat SHIP1 around the phosphorylation site of Tyr1020 (PL(p-Y)GS)
Target	Phospho-SHIP1 (Tyr1020)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	134 kDa
Uniprot ID	P97573
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung lysa...