



Product Datasheet

CD15 Antibody (orb69453)


Description	Mouse monoclonal antibody to CD15
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Immunogen	Synthesized peptide of human CD15.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	59kDa
Entrez	2526
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Immunohistochemical analysis of paraffin...



Immunofluorescence anal