

Product Datasheet

SENP2 antibody (orb694717)

Description	Rabbit polyclonal antibody to SENP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	SENP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67855
Uniprot ID	Q9HC62
NCBI	NP_067640.2
Dilution Range	IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.