


Product Datasheet

SirT1 (phospho-Ser27) antibody (orb6958)

Description	SirT1 (phospho-Ser27) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SirT1 around the phosphorylation site of Ser27 (AS(p-S)PA)
Target	Phospho-SirT1 (Ser27)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58/81 kDa
Uniprot ID	Q96EB6
Dilution Range	Flow-Cyt=1ug/Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...