


Product Datasheet

SirT1 (Phospho-Ser47) antibody (orb6959)

Description	Rabbit polyclonal to SirT1 (Phospho-Ser47). This gene encodes a member of the sirtuin family of prot...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of JK whole cell l...