

Product Datasheet

SIRT6 antibody (orb6960)

Description	SIRT6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, WB
Immunogen	KLH conjugated synthetic peptide derived from human SIRT6 (121-220/355 aa)
Target	SIRT6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q8N6T7
Expiration Date	12 months from date of receipt.