



Product Datasheet

MAP2 Antibody (orb69609)

Description	MAP2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P11137
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescent
analysis of G...



Western blot
analysis of MAP2
antibody i...