



Product Datasheet

NSFL1 cofactor p47 antibody (orb69677)

Description	NSFL1 cofactor p47 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UNZ2
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...



Immunohistochemical analysis of formalin...