



Product Datasheet

Interferon alpha-8 antibody (orb69716)

Description	Interferon alpha-8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P32881
Expiration Date	12 months from date of receipt.



Western Blot analysis
of A:Raji whole ce...



Immunohistochemical
analysis of formalin...