

Product Datasheet

TMEM40 antibody (orb697196)

Description	Rabbit polyclonal antibody to TMEM40
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	TMEM40
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25495
Uniprot ID	Q8WWA1
NCBI	NP_060776.2
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.