



Product Datasheet

CLCA1 antibody (orb69755)


Description	CLCA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	A8K7I4
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431 whole ce...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...