

Product Datasheet

Thymidine Kinase 2 antibody (orb697801)

Description	Rabbit polyclonal antibody to Thymidine Kinase 2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TK2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31005
Uniprot ID	O00142
NCBI	NP_001258864.1 , NP_001166115.1 , NP_001166114.1 , NP_004605.4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.