


Product Datasheet

interleukin 17 receptor C isoform 1 precursor antibody (orb69888)

Description	interleukin 17 receptor C isoform 1 precursor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8NAC3
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
A:Hela S3
whole...