



Product Datasheet

MYPT1 (phospho-Thr696) antibody (orb6993)


Description	MYPT1 (phospho-Thr696) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MYPT1 around the phosphorylation site of Thr696 (RS(p-T)QG)
Target	Phospho-MYPT1 (Thr696)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	O14974
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /Test
Expiration Date	12 months from date of receipt.



Western blot analysis of BRL-3A (Rat)cel...



Western blot analysis of 293T (human)cel...



Flow cytometric analysis of Hela cells u...