

Product Datasheet

SARS-CoV-2 Spike glycoprotein (C1254) antibody (orb699386)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Spike glycoprotein (C1254)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic Peptide
Target	S
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	141178
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.