

Product Datasheet

DENND4C antibody [Out of stock] (orb699423)

Description

Mouse monoclonal antibody against

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-390
Tested Applications	ICC, IP, WB
Immunogen	Full Length Human Recombinant Dennd4c protein
Target	DENND4C
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-610 was sufficient for detection of DENND4C in A549, MEF cells by Western Blot analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	213 kDa
Uniprot ID	Q5VZ89
NCBI	NM_017925
Entrez	55667
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.