

Product Datasheet

FLI1 Antibody (orb699543)

Description

Recognizes a protein of 51kDa, which is

ETS family...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

A portion of amino acids 107-261 from the human protein was used as the immunogen for the FLI1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the FLI1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the FLI1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q01543](#)

Hazard

Information

This FLI1 antibody is available for research use only.

Dilution Range

ELISA: order Ab without BSA for coating

Expiration Date

12 months from date of receipt.