

Product Datasheet

AKT1 Antibody (orb699554)

Description

Recognizes a protein of 62kDa, which is

Akt family c...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

A portion of amino acids 85-189 from the human protein was used as the immunogen for the recombinant AKT1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant AKT1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the recombinant AKT1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P31749

Hazard

Information

This recombinant AKT1 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml, ELISA: order Ab without BSA for coating

Expiration Date

12 months from date of receipt.