

## Product Datasheet

### Elastin Antibody (orb699572)

## Description

Elastin is a polymeric protein and is the main

arteries.It ...

## Species/Host

Mouse

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

ELISA, IHC-P

## Applications

## Immunogen

Recombinant full-length human Elastin protein was used as the immunogen for the Elastin antibody.

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

## Storage

Store the Elastin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

Optimal dilution of the Elastin antibody should be determined by the researcher.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

## Isotype

Mouse IgG1, kappa

## Clonality

Monoclonal

## Purity

Protein G affinity chromatography

## Uniprot ID

**P15502**

## Hazard

## Information

This Elastin antibody is available for research use only.

## Dilution Range

ELISA: order Ab without BSA for coating,Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

## Expiration Date

12 months from date of receipt.