

Product Datasheet

Cyclin D1 Antibody / CCND1 (orb699719)

Description

Recognizes a protein of 36kDa, identified as

cycle regulat...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A portion of amino acids 200-295 from the human protein was used as the immunogen for the recombinant Cyclin D1 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant Cyclin D1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the recombinant Cyclin D1 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P24385

Hazard
Information

This recombinant Cyclin D1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.