



Product Datasheet

ALPPL2 antibody (orb69977)


Description	ALPPL2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10696
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 whole ce...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...