



Product Datasheet

CLN2 antibody (orb69979)

Description	CLN2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14773
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431, B:H1299...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of SW480 xen...