

Product Datasheet

Amyloid beta Antibody / APP (orb699799)

Description

Proteolytic cleavage of the Amyloid protein

	Amyloid protein. T...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF, IHC-P
Immunogen	A portion of amino acids 578-680 from the human protein was used as the immunogen for the Amyloid Beta antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Amyloid Beta antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Amyloid Beta antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2c, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P05067
Hazard Information	This Amyloid Beta antibody is available for research use only.
Dilution Range	Immunofluorescence: 1-2ug/ml, Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.