



Product Datasheet

TNAP antibody (APC) (orb699836)

Description	Mouse monoclonal antibody against TNAP conjugated to APC
Reactivity	Canine, Goat, Human, Monkey, Porcine, Sheep
Conjugation	APC
Tested Applications	FC
Immunogen	WERI-RB-1 retinoblastoma cell line
Target	TNAP
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P05186
Entrez	249
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of TNAP positive HeLa cells (...)



Flow cytometry surface staining patterns...