


Product Datasheet


CD36 antibody Biotin) (orb699874)

Description	Mouse monoclonal antibody against CD36 conjugated to Biotin
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Platelets
Target	CD36
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P16671
Entrez	948
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of thrombocytes stained anti-

...



Flow cytometry surface staining pattern ...