


Product Datasheet


CD2 antibody (orb699888)

Description	Rat monoclonal antibody to CD2
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IP
Immunogen	murine thymocytes
Target	CD2
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Rat IgG2b lambda
Clonality	Monoclonal
Purity	Purified by protein-G affinity chromatography.
Uniprot ID	P08920
Entrez	12481
Dilution Range	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of murine CD2 positive cells

...



Flow cytometry surface staining pattern ...