



Product Datasheet

CD2 antibody (FITC) (orb699889)

Description	Rat monoclonal antibody against CD2 conjugated to FITC
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	murine thymocytes
Target	CD2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Rat IgG2b lambda
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P08920
Entrez	12481
Dilution Range	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
surface
staining
pattern ...



Separation of
murine CD2
positive cells
...