

Product Datasheet

p27 antibody (orb699901)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Western Blot (5-10 µg/ml), Immunohistochemistry, Immunoprecipitation (co-precipitates cdk4)
Isotype	IgG1
Clonality	Monoclonal
Source	The immunogen for the anti-mouse p27 was a his-tagged recombinant full length mouse p27 protein recognizing a 27kDa protein.
Uniprot ID	P46414, P46527, O08769, Q6SLL5
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.