

Product Datasheet

NGF antibody (orb699903)

Description	Mouse monoclonal antibody to NGF
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	EIA, FA, IP
Preservatives	Provided as solution in phosphate buffered saline, no preservative
Form/Appearance	Provided as solution in phosphate buffered saline, no preservative
Concentration	See vial for concentration
Storage	Store at -2°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for immunoprecipitation, EIA and for neutralizing NGF in functional studies. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG1
Clonality	Monoclonal
Source	Hybridoma produced by the fusion of splenocytes from BALB/c mice immunized with human recombinant nerve growth factor and mouse myeloma cells.
Uniprot ID	P13600, P01138, P01139
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.