

Product Datasheet

BrdU Immunohistochemistry Assay Kit (orb699905)

Description	The BrdU Cell Immunohistochemistry Kit contains all ingredients (antibodies, buffers and detection r...
Tested Applications	IF, IHC-Fr, IHC-P
Storage	Upon receipt, store entire kit at -20°C. Once the kit is thawed, you may keep it at 4°C for 5 days. For long-term storage, it is recommended you aliquot and freeze the components at -20°C, particularly the Streptavidin-HRP Conjugate, the Prediluted Biotinylated Sheep anti-BrdU Detector Antibody and the 4x Trypsin Concentrate.
Note	For research use only
Application notes	Immunohistochemistry on formalin fixed, paraffin embedded, as well as frozen tissue sections. Furthermore this kit can be used for the immunocytochemical staining of cultured cells on glass slides or cell culture plates. Optimal concentration of the primary antibody should be evaluated by serial dilutions.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	Please enquire.