


Product Datasheet

Cytokine receptor-like factor 1 antibody (orb69992)

Description	Cytokine receptor-like factor 1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	075462
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
A:HepG2
whole c...