

## Product Datasheet

### Goat anti-Rabbit IgG (H&L) antibody (FITC) (orb699967)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to FITC                                                                                                                                                                                                                                                                                                                                                          |
| <b>Species/Host</b>      | Goat                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Conjugation</b>       | FITC                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Form/Appearance</b>   | Liquid                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | Store at 2-8°C. Shelf life of this product is one year from date of receipt.                                                                                                                                                                                                                                                                                                                                                            |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Application notes</b> | <p>This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative.</p> |
| <b>Clonality</b>         | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | Antigen Immunogen affinity purified IgG, F(ab)'2 fragment.                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Dilution Range</b>    | 1 : 50-1 : 5000 (ICC) (Flow cyt), 1 : 20-1 : 2000 (IHC)                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                         |