

Product Datasheet

Rabbit anti-Guinea pig IgG (H&L) antibody (FITC) (orb700220)

Description

Rabbit polyclonal antibody against

conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Immunogen

Purified guinea pig IgG (H&L)
NP_001166520

Form/Appearance

Liquid

Storage

Store at 2-8°C. The shelf life is 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Note

For research use only

Application notes

Purity of this preparation is > 95% based on SDS-PAGE. Antibody concentration is 1.0 mg/ml (E 1% at 280 nm = 13.0). Antibody is supplied in 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2.1 % (w/v) B, Protease/IgG free. Contains 0.05% (w/v) sodium azide as preservative of bacterial growth. Based on immunoelectrophoresis, this antibody reacts with: heavy chains on guinea pig IgG, light chains on all guinea pig immunoglobulins. Based on immunoelectrophoresis, no reactivity is observed to: non-immunoglobulin guinea pig serum proteins.

Clonality

Polyclonal

Purity

Immunogen affinity purified rabbit IgG.

Dilution Range

The optimal working dilution should be determined by the investigator,
Suggested starting dilution(s): 1 : 20-1 : 2000 for most applications

Expiration Date

12 months from date of receipt.