

Product Datasheet

Goat anti-Rabbit IgG (H&L) antibody (Biotin) (orb700365)

Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to Biotin
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	IA
Immunogen	Purified Rabbit IgG (H&L)
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. No reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG, F(ab)'2 fragment.
Dilution Range	The optimal working dilution should be determined by the investigator
Expiration Date	12 months from date of receipt.