

Product Datasheet

Goat anti-Dog IgG (H&L) antibody (FITC) (orb700461)

Description	Goat polyclonal antibody against Goat anti-Dog IgG (H&L) conjugated to FITC
Species/Host	Goat
Reactivity	Canine
Conjugation	FITC
Tested Applications	FC, IL
Immunogen	Purified Dog IgG, whole molecule,
Form/Appearance	Liquid
Storage	Store at -20° C.
Note	For research use only
Application notes	Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Dog IgG (H&L). Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free. 0.05% (w/v) Sodium Azide is added as a preservative. Based on IEP, this antibody reacts with: heavy (γ) chains on dog IgG light chains on all dog immunoglobulins. Based on IEP, no reactivity is observed to: non-immunoglobulin dog serum immunoglobulins.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	1 : 20-1 : 2000 (Flow cyt), (IL)
Expiration Date	12 months from date of receipt.