

Product Datasheet

**Goat anti-Rabbit IgG (H&L) antibody [Out of stock]
(orb700564)**

Description

Goat polyclonal antibody to Goat

Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	IA
Immunogen	Purified Rabbit IgG, whole molecule
Form/Appearance	Liquid
Storage	Store at 2-8°C.
Note	For research use only
Application notes	Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Rabbit IgG. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2. 0.05% (w/v) Sodium Azide is added as a preservative. Based on IEP, this antibody reacts with: heavy (γ) chains on rabbit IgG light chains on all rabbit immunoglobulins. Based on IEP, no reactivity is observed to: non-immunoglobulin rabbit serum proteins serum proteins from bovine, goat, human, mouse, or rat IgG from bovine, goat, human, mouse or rat.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	The optimal working dilution should be determined by the investigator
Expiration Date	12 months from date of receipt.