

Product Datasheet

Goat anti-Rabbit IgG Fc antibody (FITC) [Out of stock] (orb700571)

Description

Goat polyclonal antibody against

to FITC

Species/Host

Goat

Reactivity

Rabbit

Conjugation

FITC

Tested

ELISA, FC, IF, IL

Applications
Immunogen

Purified rabbit IgG Fc,

Form/Appearance

Liquid

Storage

Store at 2-8°C. The shelf life is 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Note

For research use only

Application notes

Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Rabbit IgG. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free. 0.05% (w/v) Sodium Azide is added as a preservative. Based on IEP, this antibody reacts with: heavy (γ) chains on rabbit IgG. Based on IEP, no reactivity is observed to: non-immunoglobulin rabbit serum proteins light chains on all rabbit immunoglobulins Rabbit IgG F(ab)'₂ fragment serum proteins from bovine, horse or human

Clonality

Polyclonal

Purity

Immunogen affinity purified goat IgG.

Dilution Range

The optimal working dilution should be determined by the investigator

Expiration Date

12 months from date of receipt.