

Product Datasheet

Human TNFSF11 / RANKL / CD254 Protein (orb700632)

Description

Human TNFSF11, Mouse IgG2a Fc Tag

(HEK293). It contains AA Gly 64 ...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

TNFSF11

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-mFc(mlgG2a)

Note

For research use only

Application notes

This protein carries a mouse IgG2a Fc tag at the N-terminus. The protein has a calculated MW of 47.4 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAC51762.1

Purity

95%

Source

Human TNFSF11, Mouse IgG2a Fc Tag (orb700632) is expressed from human 293 cells (HEK293). It contains AA Gly 64 - Asp 245 (Accession # AAC51762.1).

MW

47.4 kDa

NCBI
[AAC51762.1](#)
Expiration Date

6 months from date of receipt.