


Product Datasheet

Stat2 (phospho-Tyr690) Antibody (orb7017)

Description	Stat2 (phospho-Tyr690) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P52630
Expiration Date	12 months from date of receipt.



Western blot analysis of Stat2 (phospho-Tyr690) in HeLa cells.