



Product Datasheet

SUR1 antibody (orb7031)


Description	SUR1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SUR1 (301-400/1581 aa)
Target	SUR1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	175 kDa
Uniprot ID	Q09428
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Cerebrum
(Mouse); Lan...



Western
blot: Cerebellum
(Mouse); Primary:...



Western
blot: Hippocampus
(Mouse); Primary:...