


Product Datasheet

alpha synuclein (phospho-Tyr133) antibody (orb7035)

Description	alpha synuclein (phospho-Tyr133) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SNCA around the phosphorylation site of Tyr133 (EG(p-Y)QD)
Target	phospho-Alpha synuclein (Tyr133)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	Q16143
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of
Hela cell
using
phospho-...