

Product Datasheet

SARS Spike Glycoprotein antibody (HRP) [Out of stock] (orb703857)

Description

Rabbit polyclonal antibody against

conjugated to HRP

Species/Host

Rabbit

Conjugation

HRP

Tested

ELISA, ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human SARS spike glycoprotein

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

138kDa

Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200

Expiration Date

12 months from date of receipt.