

Product Datasheet

Vinculin (phospho-Tyr822) antibody (FITC) [Out of stock] (orb703973)

Description

Rabbit polyclonal antibody against

conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human Vinculin around the phosphorylation site of Tyr822

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

124kDa

Uniprot ID
Q64727, **P85972**, **P18206**, **P12003**
Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.