



Product Datasheet

Cytochrome P450 17A1 antibody (orb704167)


Description	Cytochrome P450 17A1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Cytochrome P450 17A1
Target	Cytochrome P450 17A1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	57 kDa
Uniprot ID	P05093
Dilution Range	WB=1:500-1000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat kidney
using
Cytoc...



IHC-P
analysis of
mouse
kidney
using Cyt...



IHC-P
analysis of
Human
kidney
using Cyt...