



Product Datasheet

ACE2 antibody (orb704176)


Description	ACE2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ACE2 (180-240/805 aa)
Target	ACE2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	92 kDa
Uniprot ID	Q9BYF1
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:50, IHC-F=1:50-200, IF=1:100-200, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse kidney
using ACE...



IHC-P
analysis of
rat kidney
using ACE2 ...



IHC-P
analysis of
Human
kidney using
ACE...