



Product Datasheet

ZNF423 antibody (orb704180)


Description	ZNF423 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ZNF423 (191-290/1284 aa)
Target	ZNF423
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145 kDa
Uniprot ID	Q2M1K9
Dilution Range	WB=1:500-2000, IHC-P=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat testis using ZNF423...



IHC-P analysis of mouse testis using ZNF...



Western blot: Lane 1: Cerebrum (Mouse); Lan...